



Graphic Display



Keypad 220g/0.0001g
Analytical Balance
KAB2L



SALIENT FEATURES

- ✦ Graphic Display
- ✦ Internal Calibration
- ✦ GLP / ISO Record of last stored calibration data
- ✦ EMFC Based, Transportation lock
- ✦ Memory 99 receipts can be stored
- ✦ User selectable environmental filter
- ✦ Piece Counting, Density, Check Weight, Percentage, Animal Weighing, Formulation,
- ✦ Hi-Lo limit Alarm, Date & time
- ✦ Underhook Weighing, RS232C interface
- ✦ Plastic protection cover for Display
- ✦ Multi language Option

FORMULATION:

 03-04-11 10:54:57

 Manual
 1. 31.05 g
 2. 100.02 g
 3. 26.89 g

 T = 157.96 g

 NOTE: To transmit the print of total of weights, press and keep pressed the PRINT button

Value under threshold	Value inside thresholds	Value over threshold
03-04-11 11:02:19	03-04-11 11:01:50	03-04-11 11:01:50
-----	-----	-----
Lim.1: 10.00 g	Lim.1: 10.00 g	Lim.1: 10.00 g
Lim.2: 100.00 g	Lim.2: 100.00 g	Lim.2: 100.00 g
Weight: -0.01 g	Weight: 31.08 g	Weight: 131.10 g
TEST: KO! -	TEST: OK!	TEST: KO! +++

PERCENTUAL WEIGHT

 03-04-11 11:58:39

 Perc. 100.0 %
 Weight: 18.69 g
 Refer.: 18.69 g

Analytical Balance

HSN Code - 90160010

SPECIFICATIONS

MODEL	KAB2L
Make	WENSAR
Max. (Capacity)	220g
Readability (d)	0.0001g (0.1mg)
Repeatability	0.0001g
Linearity (±)	0.3mg
Min. Capacity	10mg
Verification Scale Intervals	1mg
Calibration	Internal Calibration
Class	Class I
Display	Graphic Display with Back Light
Printout of Results	Date, Time, Instrument ID, User ID, Project ID, Weight, Signature
Response Time	≤ 4 Sec
Pan Size (mm)	80mm Ø
Tare Range	Full
Operating Temperature	15°C to 30°C
Sensitivity Drift (±)	2ppm/°C
Density	Function Available
Key Sound	Yes
Data Output Mode	Command, Continuous, Stable
Unit	g, mg, ct, oz
Baud Rate Options	2400, 4800, 9600
Interface	RS232C
Power supply	DC Adaptor, Using AC 220V / 50-60Hz
Packing Dimension (L x W x H)	515 x 400 x 525mm
Gross Weight	8kg

* Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Weighing PAN, Adaptor, Instruction Manual, Test Certificate.