














220g/0.1mg


Touch Screen **Analytical Balance**
TAB1T

SALIENT FEATURES

-  Electro Magnetic Force Compensation
-  Touch Screen display
-  Percentage Weighing
-  Multiple Weighing Units
-  Density Determination Function
-  RS232C Interface
-  Piece Counting Function
-  Underhook Weighing
-  Calibration with Internal weight
-  Time
-  Date



E - Catalogue



Touch Screen Analytical Balance

HSN Code - 90160020

SPECIFICATIONS

MODEL	TAB1T
Make	WENSAR
Max. (Capacity)	220g
Readability (d)	0.1mg
Repeatability	0.2mg
Linearity (±)	0.4mg
Min. Capacity	10mg
Verification Scale Intervals (e)	1mg
Calibration	Internal Calibration
Class	Class I
Display	5 inch Touch Screen
Response Time	3 - 5 Sec
Pan Size (mm)	80 mm Ø
Tare Range	Full
Operating Temperature	15°C to 30°C
Sensitivity Drift (±)	2ppm / °C
Density	Function Available
Key Sound	Yes
Weighing Chamber	240 x 190 x 265 mm
Unit	g, mg, ct, oz
Baud rate Options	1200, 2400, 4800, 9600
Interface	RS232C
Power Supply	AC 220V / 50-60 Hz
Packing Dimension (LxWxH)	510 x 370 x 515mm
Gross Weight	10kg

* Specifications are subject to change without notice.

* Density Determination Kit for Liquid or Solid Samples separately (Optional)

STANDARD ACCESSORIES : Main Instrument, Weighing PAN, Power Cable, Dust Cover, Instruction Manual, Test Certificate.