

Antivibration Pad



AVP304SS - Construction SS304



AVP100 - Construction MS

SALIENT FEATURES

- ✦ Antivibration Pad for Micro/ Semi Micro / Analytical Balance
- ✦ Easy to Clean.
- ✦ Good shock Absorption effect.
- ✦ No Noise.

Antivibration Pad offers a solid & stable surface developed specifically to help decrease vibration during weighing measurement. When weighing balance used with Anti Vibration pad, its allow the balance to perform with marked precision, despite air currents or movements that may cause readings to fluctuate.



E - Catalogue

HSN : 84798200

SPECIFICATIONS

MODEL	AVP304SS	AVP100
Construction	SS304	MS
Length	470mm	
Width	330mm	
Height	75mm	
Top Colour	Black	
Noise Level	Nil	
Thickness	1mm	
Construction Material	Galvanized Steel Plate	
Vibration	In Built Shock Absorption Device	
Level	Manually	
Power Supply	NA	
Packing Dimension (LxWxH)	480 x 350 x 90mm	
Gross Weight	12kg	

* Specifications are subject to change without Notice

Note : Its suitable for any make Micro, Semi Micro / Analytical Balance