



SALIENT FEATURES

- ✦ Aluminium Alloy Case
- ✦ White Back Light LCD display
- ✦ Piece Counting
- ✦ Overload alarm Level Indicator
- ✦ Stainless Steel Square PAN
- ✦ RS232C Interface

3000g/0.01g

Precision Balance

PGB3001



E - Catalogue

HSN Code - 90160010

SPECIFICATIONS

MODEL	PGB3001
Make	WENSAR
Capacity	3000g
Readability	0.01g (10mg)
Repeatability	0.01g
Linearity (±)	0.03g
Pan Size	180 x 180mm
Calibration	External
Display	LCD
Unit	g, ct, oz, gsm
Operating Temperature	15°C to 35°C
Power supply	DC Adaptor, Using AC 220V / 50Hz
Packing Dimension (L x W x H)	400 x 265 x 165mm
Gross Weight	3kg

* Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Weighing Pan, Adaptor, Instruction Manual, Test Certificate.