



**SALIENT FEATURES**

- ✦ External Calibration
- ✦ Easy to read LCD Display with backlight
- ✦ Multiple Weighing Units
- ✦ Auto Zero tracking, Low Battery Indicator
- ✦ Overload Protection Design
- ✦ S.S. Weighing Pan, Level indicator, Adjustable Fit
- ✦ High Resolution, Quick Weighing, Accurate Result
- ✦ PGB3010 with RS232C Interface

**3000g/0.01g    6000g/0.1g**  
**Precision Balance**  
**PGB3010 & 6000**



HSN Code - 90160010

HSN Code - 84238900

**SPECIFICATIONS**

MODEL	PGB3010	PGB6000
Capacity	3000g	6000g
Readability	0.01g (10mg)	0.1g (100mg)
Repeatability	0.01g	0.1g
Linearity (±)	0.03g	0.2g
RS232C Interface	Available	NA
Make	WENSAR	
Pan Size	180 x 195 mm	
Display	LCD	
Calibration	External Calibration	
Tare Range	Full	
Operating Temperature	10°C to 40°C	
Power supply	DC Adaptor, Using AC 220V / 50Hz	
Packing Dimension (LxWxH)	380 x 265 x 130mm	
Gross Weight	3kg	

\* Specifications are subject to change without notice.

**STANDARD ACCESSORIES :** Main Instrument, Pan, Adapter, Instruction Manual, Test Certificate.