



SALIENT FEATURES

- ✦ Internal Calibration
- ✦ Touch Screen display (7" full Color HD)
- ✦ Stainless steel large size square PAN
- ✦ Height adjustable feet
- ✦ Date & Time Setting
- ✦ Piece Counting Function, Accumulation
- ✦ Check Weighing, Under Weighing, %, Density Function
Net gross weight / sum
- ✦ Overload Alarm, Connect PC and Printer

30kg/0.1g

TableTop Balance

TTB30M



E - Catalogue

HSN Code - 84238900

SPECIFICATIONS

MODEL	TTB30M
Capacity	30kg
Readability	0.1g
Repeatability	0.1g
Linearity (\pm)	0.3g
Make	WENSAR
Pan Size	305 x 230mm
Calibration	Internal Calibration
Display	Touch Screen
Operating Temperature	5°C to 35°C
Unit	g, oz, ct...
Application Mode	Piece Counting, Check Weighing, Percentage, Density
Maximum Consumed Power	5W
Interface	RS232C / USB
Power supply	DC Adaptor, Using AC 220V / 50Hz
Packing Dimension (L x W x H)	465 x 370 x 210mm
Gross Weight	7kg

* Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Weighing PAN, Power Adaptor, Instruction Manual, Test Certificate.