



Overhead Stirrer
HOS20L

SALIENT FEATURES

- ✦ Digital LCD Display.
- ✦ 100 ~ 2500rpm speed.
- ✦ Maximum Volume 20L.
- ✦ Stirring has three modes - Clockwise,
- ✦ Anticlock wise & Alternating.
- ✦ Stainless Steel Head.
- ✦ Digital type bushless motor.



E - Catalogue

HSN : 84798200

SPECIFICATIONS

MODEL	HOS20L
Max. Volume	20L
Display	LCD
Speed(r.p.m)	100 ~ 2500
Working Time	Continuous
Structure Exterior	Plastic shell
Motor power	60W
Setting & Speed mode	Digital
Timer	Yes
Max Viscosity	10000 mpas
Plate size	200 × 300 mm
Plate load	50kg
Power Supply	Using Adaptor, AC 220V / 50-60Hz
Packing Dimension & Gross Weight	400 × 265 × 185mm & 8kg

* Specifications are subject to change without Notice

Standard Accessories : Main Instrument, Power Cable with adaptor, Clamp, Washer, Chuck Key, Mixing paddle For Sample, Stand with Extension , Instrument Holding shaft, Spanner Base Plate, Instruction Manual, Test Certificate .