



SALIENT FEATURES

- ✦ 20L Maximum Stirring Volume (1L to 20L)
- ✦ LCD Display, PT1000 Temperature Sensor
- ✦ Speed & Timing Adjustable
- ✦ Bushless Motor Type, Low Noise
- ✦ Function of forward and reverse rotation.
- ✦ High Temperature Resistance
- ✦ Top plate, Cast Aluminium with Ceramic Coating
- ✦ Automatically memorise the set working parameters
- ✦ Aluminium alloy base

Ceramic Hotplate Magnetic Stirrer

LMMS20LC



E - Catalogue

HSN Code : 84798200

SPECIFICATIONS

MODEL	LMMS20LC
Make	LABMAN
Stirring Volume (H ₂ O)	5L Standard (1L -20L)
Speed range	100-2000RPM (Adjustable)
Plate Size	170 x 170 mm
Top Plate Type	Cast Aluminium with Ceramic Coating
Speed Control Accuracy	±10 RPM
Display	LCD
Temperature Sensor	PT1000
Time Setting	Yes
Max. Temperature on Plate	350°C
Power Wattage	950 W
Power supply	AC 220V / 50Hz
Packing Dimension (LxWxH)	360 x 295 x 170mm
Gross Weight	4kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, ATC Probe, Teflon Paddle, Clamp, Bracket, Power cord, Instruction Manual, Test Certificate.