



Multiposition Magnetic Stirrer

MMS15S - Stirring Facility

MMS12H - Stirring & Heating Facility

- ◆ LCD display, 4 Position Magnetic Stirrer, Max capacity 500ml for each Position.
- ◆ Timing and locking function, Chemical & corrosive resistance
- ◆ Bushless motor, Mixing Direction in Clockwise

Multiposition magnetic stirrers are ideal for laboratories that require simultaneous stirring of multiple samples.



E - Catalogue

HSN Code : 84798200

SPECIFICATIONS

MODEL	MMS15S	MMS12H
Speed	50 - 1500 RPM	50 - 1200 RPM
Plate Material	PET Anti-Corrosion Plate	304 Stainless Steel Plate
Temperature Range	NA	RT+5°C to 120°C
Make	LABMAN	
Control	Asynchronous Control (Synchronization can be set)	
Number of Stirring Points	4	
Max Stirring Volume	4x500ml	
Plate Size	100 x 110mm (each)	
Teflon Paddle Size	30mm	
Display	LCD	
Time Setting	Yes	
Power Supply	AC 220W / 50-60Hz	
Ambient Temperature	5 ~ 40°C	
Packing Dimension (LxWxH)	380 x 290 x 115mm	380 x 300 x 115mm
Gross Weight	2kg	3kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Magnetic Stirring Bar, Power Cord, Instruction Manual, Test Certificate.