




**SALIENT FEATURES**

- ◆ Stainless Steel Inner Tank
- ◆ Built-in heating tube - safe experiment
- ◆ Transparent long cover for better view inside
- ◆ PID temperature control
- ◆ Bushless Motor
- ◆ Anti - Dry Burning Function

**Waterbath Magnetic Stirrer**

**WBMS3**



E - Catalogue

HSN Code : 84798200

**SPECIFICATIONS**

MODEL	WBMS3
Make	LABMAN
Maximum Mixing Volume	3L (2L flask or 3L Beaker)
Display	LCD
Speed	100RPM~1200 RPM
Increment	50 RPM
Tank Size	200 x 200 x 130mm (depth)
Temperature Range	RT +5°C~100°C
Temperature Resolution & Stability	0.1°C & ±1°C
Control Mode	Digital Control Panel Key
Time Setting	Yes
Power	800W
Power supply	AC 220V / 50-60Hz
Packing Dimension	300 x 300 x 250mm
Gross Weight	4kg

★ Specifications are subject to change without notice.

**STANDARD ACCESSORIES:** Main Instrument, ABS Top Cover, Stirring Bar, Power Cord, Instruction Manual, Test Certificate.