



Column Oven LMCO1

SALIENT FEATURES

- Large White LCD Display, Clear Text View
- Chamber is equipped with a three slot double column base
- Analytical column and protective column is within 340mm
- Factory set temperature at 40 °C
- Self Tuning Function (PID) and real time Indication
- Temperature Setting & Control Operation
- Integral thickened alloy, which is stable and reliable




HSN Code : 84178090

SPECIFICATIONS

MODEL	LMCO1
Make	LABMAN
Temperature range	RT+5°C~80°C
Constant Temperature accuracy	0.1°C
Constant Temp. chamber size(LxWxH)	80 x 38 x 400mm
Length of Column	340mm
Power Wattage	120W
Power supply	AC 220V / 50Hz
Packing Dimension (LxWxH)	550 x 210 x 215mm
Gross Weight	3kg

★ Specifications are subject to change without notice.



Real-time Temperature
Setting Temperature
Set/confirm Button

Shift Key
Down Key
Up Key

STANDARD ACCESSORIES : Main Instrument, Power Adaptor, Instruction Manual, Test Certificate.