



Orbital Shaker

OS20

SALIENT FEATURES

- ✦ Digital Display
- ✦ Speed Control knob to adjust the RPM.
- ✦ Speed range is 40-200 rpm.
- ✦ Timing knob can set the duration of operation.
- ✦ Start/Pause Switch for start or pause the operation.
- ✦ Universal Tray Platform for holding various containers.
- ✦ Bushless motor.
- ✦ Suitable for Various Containers

Orbital shakers are used for mixing, aerating, and incubating samples in laboratory, particularly in cell culture, biochemical reactions, and sample preparation.



E - Catalogue

HSN Code : 84798200

SPECIFICATIONS

MODEL	OS20
Make	LABMAN
Display	LCD
Speed Range	40-200rpm
Plate Type	Universal Tray Platform
Plate Size	250x250x50 mm
Orbital Motion	Circular
Max Load	2kg
Orbit diameter	20mm
Motor Type	Bushless Motor
Timing Range	Yes
Ambient Temperature	5°C to 40°C
Relative Humidity	80%
Power Supply	AC 220V/50-60Hz
Packing Dimension (LxWxH)	410 x 300 x 170mm
Gross Weight	5kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Power cord, O Ring, Instruction Manual, Test Certificate