



pH Meter - Economic
LMPH9

SALIENT FEATURES

- ◆ Automatic Temperature Compensation
- ◆ There is a microprocessor chip inside the meter and the appearance of meter looks nice with user friendly operation. Some of the features are as under
- ◆ Large LCD display with blue background simultaneous display of pH and temperature
- ◆ It has intelligent functions for calibration, data storage / Memory (upto 25 prints)
- ◆ Automatically display the percentage of slope after calibration is over.
- ◆ Equipped with electrode holder and calibration buffer solution, which is more convenient to use.



E - Catalogue

HSN Code : 90275090

SPECIFICATIONS

MODEL	LMPH9
Make	LABMAN
Range	0 to 14.00 pH
pH Resolution	0.1 pH / 0.01 pH
Accuracy	± 0.01 pH
Calibration Points	up to 3 points
pH Buffer options	USA (pH 4.01, 7.00, 10.01)
mV Range	±1999 mV
Resolution & Accuracy	± 1 % FS ± 1 digit & 1mV
Temperature Range	0°C to 100°C
Accuracy	± 0.5°C
Resolution	0.1°C
Temperature Compensation	Automatic
Connector	BNC
Display	LCD
Memory	Upto 25 prints
Power Requirements	DC Adaptor, Using AC 220V,50/60Hz
Packing Dimensions (Lx W x H)	300 × 250 × 130mm
Weight	1kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : pH Meter, pH Electrode, pH Buffer - 4.01, 7.00, Power Adapter, Electrode Stand, ATC probe, Instruction Manual, Test Certificate.