



## Vortex Mixer

### LMVM60D

### SALIENT FEATURES

- ◆ LCD Display with white back light
- ◆ Speed Timer (1~99h 59m 59s)
- ◆ Operating Mode - 2 Types - Auto run sensor, continuous run used for various Mixing application
- ◆ Supply with Standard Head
- ◆ Soft silicone foot pads
- ◆ Shockproof, Brushless Motor
- ◆ Aluminium Alloy Base stable body

Digital display vortex mixer is powered by a brushless motor and features an eccentric structural design.

This innovative mixer effectively blends and oscillates liquids within a container. It finds extensive applications in various fields, including biochemistry, genetic engineering, and medicine for conducting experiments



E - Catalogue

HSN Code : 84798200

### SPECIFICATIONS

MODEL	LMVM60D
Make	LABMAN
Speed	100-2800rpm (adjustable)
Display	LCD with white back light
Display Function	Speed & Time
Adjustment method	Keypad
Operation Mode	2 Types - Auto run sensor, continuous run
Oscillation mode	Circumference
Base material	Aluminium alloy
Shaking Movement	Orbital
Orbital diameter	4mm
Ambient Temperature	5°C-35°C
Humidity	≤ 80%
Input power	60W
Power supply	AC 220V,50/60Hz
Packing Dimension (LxWxH) & Gross weight	240 x 200 x 190mm & 4kg

★ Specifications are subject to change without notice.

**STANDARD ACCESSORIES:** Main Instrument with Standard Head, Power Cable with Adaptor, Instruction Manual, Test Certificate.