



Multi Tube Vortex Mixer
MTVM50

SALIENT FEATURES

- ◆ 50 samples can mix simultaneously
- ◆ LED Display to show both timing and speed.
- ◆ Timer Function, Brushless Motor
- ◆ Horizontal Orbital Shaking
- ◆ 5 Types of Sponge Test tube racks
- ◆ Accurate Speed Control
- ◆ Impulse & Continuous mode

A multi-tube vortex mixer is a lab device that simultaneously mixes multiple tubes, ensuring efficient and uniform blending of samples. It's ideal for high-throughput laboratory applications, saving time and increasing productivity.



E - Catalogue

HSN Code :84798200

SPECIFICATIONS

MODEL	MTVM50
Make	LABMAN
Maximum Load Capacity	4kg
Operation Mode	2 Modes
Display	LED
Speed Range (rpm)	500-3000rpm
Replaceable Sponge Rack	5 types
Test tube Hole Size	6mm, 8mm, 15mm, 25mm, 29mm
Maximum No. of Hole	50
Speed Adjustment Increment (rpm)	1rpm
Amplitude	4mm
Timer	Yes
Time Display Accuracy	1Sec
Time Adjustment Increment	1min
Power (W)	60W
Power Supply	AC 220V / 50-60Hz
Instrument Dimension (LxWxH)	385x 235x 420mm
Packing Dimension (LxWxH)	450x 295x 495mm
Gross Weight	17kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Powercord, Instruction Manual, Test Certificate.

STANDARD ACCESSORIES



6mm x 50 Holes



15mm x 50 Holes



8mm x 50 Holes



25mm x 15 Holes

Sponge Test Tube Rack

- 6mm x 50 Holes
- 8mm x 50 Holes
- 15mm x 50 Holes
- 25mm x 15 Holes
- 29mm x 15 Holes



29mm x 15 Holes