



**SALIENT FEATURES**

- ◆ TRM6 is efficient & gentle mixing in various laboratory
- ◆ Digital key control, LCD Display
- ◆ Max Capacity: 10 x 50ml Centrifuge tube/0.5kg
- ◆ Adjustable speed, Bushless Motor,
- ◆ Silicone rollers. No. of. Roller-6,
- ◆ Roller Length-280mm, Antilock Design



A tube roller mixer is a lab device that gently rotates tubes to mix and blend samples, ensuring uniform distribution without foaming or damage.

**Tube Roller Mixer  
TRM6**

**Applications:** It is used in Scientific Research Institutions, Medical Schools, health care and other units, Molecular Biology, Virology, Microbiology, Pathology, Immunology and Laboratories.



E - Catalogue

HSN Code : 84798200

**SPECIFICATIONS**

MODEL	TRM6
Make	LABMAN
Speed	20-80rpm
Ambient Temperature	+5°C - +40°C
Relative Humidity Range	≤80%
Type	Rotating Shaking
Max. Capacity	10x50ml Centrifuge tube/0.5kg
Control Mode	Digital key
Display	LCD
Time Setting	Yes
Motor	Bushless
Rollers	6 Rollers
Roller Size	φ6-34mm
Roller Length	280mm
Power Supply	AC 220V / 50 - 60Hz
Packing Dimension (LxWxH)	510 x 305 x 145mm
Gross Weight	3kg

★ Specifications are subject to change without notice.

**STANDARD ACCESSORIES :** Main Instrument, Power Cable , Power cord, Instruction Manual , Test Certificate